

## HYDROMETER TEST SET

### ASTM D422 / AASHTO T88

This equipment is used to determine the quantitative distribution of very fine particle size in soils such as clay and silt. Set includes:

- 1 piece of Transparent Water Bath
- 1 piece of 151 H Hydrometer
- 1 piece of 151 H Hydrometer
- 1 piece of 0 - 500C Glass Thermometer
- 1 piece of 250 ml Beaker
- 1 piece of Glass Baget
- 1 piece of High Power Stirrer
- 6 pieces of 1000 ml Capacity Hydrometer Measures



#### SPECIFICATIONS:

Product Code : BSO-205  
 Weight (Set) : 19 kg  
 Dimensions : 300 x 630 x 530 mm

## END OVER END SHAKER

### BS 1377-2

It is used to determine the specific gravity of soils. It rotates two glass jars at approx. 50 rpm. Glass jar is 1 litre in capacity with rubber bung.

#### SPECIFICATIONS:

Product Code : BSO-204  
 Power : 220 V,50 Hz  
 Weight : 13 kg  
 Dimensions : 400 x 870 x 330 mm

